

One-Step SmartPee® Cup

» One-Step to on-site drug testing



AMP AMP AMP AMP	Amphetamine
BAR BAR BAR BAR	Barbiturate
BZD BZD BZD BZD	Benzodiazepine
COC COC COC COC	Cocaine
MET MET MET MET	Methamphetamine
MTD MTD MTD MTD	Methadone
OPI OPI OPI OPI	Opiates
OXY OXY OXY OXY	Oxycodone
TCA TCA TCA TCA	Tricyclic Antidepressants
PCP PCP PCP PCP	Phencyclidine (PCP)
THC THC THC THC	Marijuana
MDMA MDMA MDMA MDMA	Ecstasy
PPX PPX PPX PPX	Propoxyphene
BUP BUP BUP BUP	Buprenorphine

» **PRODUCTS** that care
for you «

- 15-panel test
- Available with adulterant and/or alcohol tests
- 8 oz. cup
- 9% larger opening
- Unbreakable Cup surface
- Results in 5 minutes
- 15 individual chambers eliminate drug interpretation confusion
- Tamper-proof strip insert
- Leak-proof transport
- No assembly, no card dipping, no moving parts
- No tilting to activate
- No cross contamination

Preparation, Specimen Collection, Test Procedure and Interpretation of Results



Collectors should be trained on the procedures for operating the tests. Do not break the seal of the protective pouch until ready to begin testing. Visually inspect the foil package to ensure it has not been compromised before beginning the test. If the package is not intact, the integrity of the test may be compromised.

1. Fresh urine specimens should be collected directly into the collection cup and do not require any special handling or pre-treatment.
2. The SmartPee® Test Cup employs a temperature strip to validate that the urine is a freshly voided specimen. This device should be checked immediately after collection.

Note: Urine specimens can be transferred from a urine collection container into the SmartPee® Test Cup, if necessary.

3. Tear open the foil pouch and remove the test cup.
4. Issue the device to the individual to be tested. Instruct the donor not to flush the toilet or turn the water on and to open the door when they are done.
5. Have the individual urinate directly into the SmartPee® Test Cup. The specimen should reach the level of the 30 ml line on the test label. The minimum specimen volume is 10 ml!
6. Check temperature (normal range between 90 °F - 100 °F or 32 °C - 38 °C) and observe specimen for foreign material and/or discoloration.
7. Secure container by tightening lid as far as it will turn, apply security seal and have donor/professional both sign and date seal.

Drug test results

PRESUMPTIVE POSITIVE AND CONFIRM

One colored line only is visible in each Control Zone. No colored line whatsoever appears in the appropriate Test Zone, indicating a presumptive positive result for the corresponding drug of that specific Test Zone. Send this urine specimen to a certified laboratory for confirmation.

NEGATIVE

Two colored lines are visible in each Control Zone and in the appropriate Test Zone, indicating that the concentration of the drug is below the detection limit of the test or the drug is not present.


RETEST

No lines: If a colored line is not visible in the Control Zone, the test is invalid. Another test should be run to re-evaluate the specimen. Each strip in the SmartPee® Test Cup is read and functions independently. An invalid result on one test strip does not invalidate other results derived from the same device.

Note: There is no meaning attributed to line color intensity or width. Any evidence of a line should be considered a line. In order to prevent any incorrect results, the test results should not be interpreted after 10 minutes.

If confirmation is required, secure specimen in appropriate shipping container and send to Confirmation Lab.

bluepoint® MEDICAL sales partner



SmartPee® is a registered trademark of bluepoint® MEDICAL.

For further information, please visit www.bluepoint-medical.com

bluepoint MEDICAL



Just collect specimen right into the cup.



Peel the privacy label.



Read results in just 5 minutes.

